

# Medical Consumables

- Tube (Infiltration, Suction, HVP, Microaire)
- Skin Stapler
- Cautery Pad, Cautery Pencil
- Dressing



## Atofuse™

### Infiltration Tubing

This especially soft formulated tubing will work in most infiltration pumps. 9' operational length with a silicone tubing insert to increase flow, eliminate tubing creep and make changing tubing more convenient. Complete with luer lock and bag spike fittings. Packaged sterile, latex-free. 10/Box.



REF AFTU-DS



## Suction tube

### Aspiration Tubing

Hi-Vac™ Aspiration Tubing is the long time market leader because of its superior quality. Hi-Vac™ tubing is highly resistant to kinking or collapsing. Available in 8' lengths for standard aspirators, or 12' lengths for PSI-TEC™ systems. Packaged sterile, latex-free. 10/Box.

REF 1240100

## Suction Tub Micro Aire Tube

The MicroAire Suction Tube, 3 meters in length, is specifically designed to integrate seamlessly with the Power-Assisted Liposuction (PAL) system. Engineered for optimal performance, it delivers precise suction capabilities, enhancing the efficiency and reliability of the PAL device across various clinical applications.



REF ML1240120

HVP™ Re-Injection Tubing,  
Sterile (10pcs)

## HVPM RE-INJECTION TUBING SET

Proprietary HVPTM Re-Injection Tubing is 12ft in length and specifically engineered for use with the HVP™ Auto-Graft Infuser. The HVPTM Re-Injection Tubing is required for the correct operation and functionality of the HVP™ system features. This tubing is precision-extruded and optically inspected to insure top performance with the pump head. Other commercial tubing may not meet, or perform, to specifications. Inspect the tubing package prior to opening for any damage and discard entire packaging and contents if the packaging appears compromised.



REF 1240160

## Disposable Skin Stapler

Features:

1. single use, avoid cross infection.
2. Easy use, save time
3. Minimal trauma, suitable for different skin incision suture.
4. Special staple remove, reduce pain.

Application : Used for trauma and incision suture of surface layer skin during the surgery operation.



REF HASPF-35W



## Cautery pencil

Advanced high-precision cautery pencil for reliable performance.  
Provides controlled coagulation with minimal bleeding and tissue damage.  
Suitable for surgical procedures, dermatology, and aesthetic clinics

REF

## Disposable Adult Bipolar Cautery Pad

A single-use bipolar cautery pad designed for safe and reliable electrosurgical procedures. It provides secure adhesion, even current distribution, and high patient safety while reducing the risk of cross-contamination.

Applications:

- \* General surgery
- \* Electrosurgical procedures
- \* Hospital operating rooms and surgical centers



REF HASPF-35W